







**COMPATIBILITY LIST**

Brand	Model	Series									Connector
		F7950	F7950W	F7750	F7750W	F7950P(1)	F7950PW(1)	F7950PL	F7950WPL	F7750P(1)	
PHILIPS	Heartstream ForeRunner AED (E, S, EM)	●	●	●	●	-	-	-			  F7950 F7950W F7750    F7950P F7950PL F7750P F7950PW   
	HeartStart FR AED (940010XX, 940020XX, 940030XX)	●	●	●	●	-	-	-			
	Heartstream/HeartStart FR2 AED (M3840A, M3841A), FR2+ AED (M3860A, M3861A)	●	●	●	●	-	-	●	●	-	
	HeartStart FR3 AED (ECG 861389, Text 861388) *	●	-	●	-	-	-		-	-	
	CodeMaster 100 (M2475B), XE (M1724A), XL (M1723A/B), XL+ (M1722A/B) **	●	●	●	●	●	●	-	-	●	
	HeartStart 4000 (M5500B) ***	●	●	●	●	●	●	-	-	●	
	Heartstream/HeartStart XL (M4735A), XLT (M3500B) ***	●	●	●	●	●	●	-	-	●	
	HeartStart MRx (M3535A, M3536A) ***	●	●	●	●	●	●	-	-	●	
	HeartStart XL+ (861290) ***	●	●	●	●	●	●	-	-	●	
	Efficia DFM100 (866199) ***	●	●	●	●	●	●	-	-	●	

F7950PL, F7950WPL: Reduced energy: defibrillator output of 150J is reduced to 50J

\* Electrodes will operate with device for therapy, but cannot be left attached while FR3 is in standby mode or while turned off due to self-test activity of the device

\*\* With hands-free pad cables M1750A/B or M2472A and adapter 05-10200

\*\*\*With hands-free pad cables M3508A